

QuickVet[®]
Analyzer

YOUR PLATFORM FOR BETTER DIAGNOSTICS




MICRO FLUIDIC TECHNOLOGY

ADD SAMPLE . READ RESULT

zoetis

THE PLATFORM FOR THE FUTURE “IN CLINIC” TESTING



The QuickVet® Analyzer is a novel veterinary diagnostic testing platform intended for in vitro testing of a range of different indications in companion animals.

The measuring principle used in the QuickVet® Analyzer is applicable to testing for a large range of indications within infectious diseases, parasites and parasite transmitted diseases, toxicology, hormones, oncology, hemostasis and cardiology in whole blood or plasma. A number of new test cartridges are under development and will be introduced over the next years, giving you the full benefits of a single testing platform.

A SYSTEM FOR BETTER DIAGNOSTICS BASED ON BRAND NEW TECHNOLOGY

The QuickVet® Diagnostic System consists of a tabletop analyzer and single-use disposable test cartridges based on Capillary Driven Micro Fluidic technology.

It is designed to perform tests using whole blood or plasma samples in veterinary clinics at the patient side in just a few minutes. Operating the instrument requires minimum training, just insert a test cartridge, add sample and read the result.

QUICKVET® ANALYZER

Built on state-of-the-art lab-on-a-chip technology the QuickVet® Analyzer features:

- Fast and accurate patient side testing
- Simple and user-friendly operation
- Very small sample volume requirement
- No or limited sample preparation
- No regular service necessary
- Instrument with small footprint
- Large 7” color touch screen
- Up to 1000 test results can be stored and retrieved
- Allows for Patient and Sample ID input
- Continuous internal self-test during operation.



QuickVet® is a registered trademark of Zoetis Denmark Aps.